Thopaz+™

Cardiothoracic Drainage

EN Instructions for use

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SV Bruksanvisning

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
Thopaz+ University

At the web address below you will find short how-to-use videos for Thopaz+.

www.medela.com

Thopaz+ University is designed to support you. Each video lasts no longer than 60 seconds and includes text bullet points and voice commentary.

Thopaz+ University can be found at www.medela.com via the "University" tab. It is available on all formats including desktop, smartphone and tablet. We recommend saving Thopaz+ University as a shortcut on desktop, smartphone and tablets for even quicker access.

Table of Contents

Congratulations	6
1 Warnings and safety instructions	7
Safety instructions.....	10
2 Power supply and battery operation	11
Backup battery.....	11
Charge battery.....	11
Battery low.....	12
Battery empty.....	12
3 Description	13
Introduction.....	13
Intended use/indications.....	13
Contraindications.....	13
Intended user.....	13
Intended patient population.....	13
Important note.....	13
4 Overview	14
Main elements of Thopaz+.....	14
Display.....	15
How to activate a button.....	15
Tubings.....	16
Canisters.....	18
5 Additional information	19
Definition of vacuum.....	19
Regulated pressure.....	19
Safety pressure limit.....	19
Tube flushing.....	19
Safety chamber in canister.....	19
Fluid measurement.....	20
Limitations.....	20
Overflow protection/bacteria filter.....	21
Light sensor.....	21
6 Installation	22
Check initial delivery.....	22
Initial startup.....	23
7 Preparation for use/operating instructions	25
Checks before use.....	25
Connect tubing.....	25
Snap in canister.....	26
Using 2l canister.....	26
Switch Thopaz+ on.....	27
Carry out functional check.....	28
Check therapy progress.....	29

Main display.....	29
Reset fluid display.....	30
History.....	31
Catheter check.....	33
Switching Thopaz+ off.....	34
Transferring data to PC with ThopEasy+.....	35
8 Replace canister.....	36
Canister with solidifier: activation of the solidification process.....	37
9 Changing pressure during operation.....	38
Physiological pressure.....	38
10 Changing settings.....	39
Factory settings.....	39
Change settings.....	40
Save a changed parameter as the new standard.....	40
11 Taking a drainage sample.....	41
12 Troubleshooting and alarms.....	42
Description.....	42
13 General cleaning guidelines.....	45
General notes.....	45
Disposables.....	45
14 Warranty, maintenance and checks.....	47
Warranty.....	47
Maintenance.....	47
Routine check.....	47
Instructions for FPT test.....	47
Safety check.....	48
15 Disposal.....	49
16 Accessories.....	49
17 Technical specifications.....	50
18 Signs and symbols.....	51
19 Technical documentation.....	413
Appendix A.....	418
Appendix B.....	420

These instructions for use belong to firmware 1.01.

Congratulations

With Thopaz+, you have acquired an innovative digital chest drainage system that sets new standards in therapy and chest drainage management.

Thopaz+ has an electronic measuring and monitoring system with optical and acoustic status indications. The device is a dry system, which means no fluids are necessary for operation. Important information regarding the course of therapy is indicated digitally and graphically on the display. Upon completion of therapy, it can be transferred to a PC.

The handy size, the extreme quietness of operation and the light sensor regulating the brightness of the screen are very pleasant for the patient as well as those in the vicinity.

1 Warnings and safety instructions

WARNINGS

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTIONS

Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

Safety related tip

Indicating useful information about the safe use of the device.

Thopaz+ is approved exclusively for the use as described in these instructions for use. Medela can only guarantee the safe functioning of the system when Thopaz+ is used in combination with the original Thopaz+ accessories (canisters, tubings, carrying strap, mains adapter, docking station – see chapter 16 – Accessories/Appendix A).

Additional equipment connected to medical electrical equipment must comply with the respective IEC or ISO standards (e.g. IEC 60950 for data processing equipment). Furthermore all configurations shall comply with the requirements for medical electrical systems (see IEC 60601-1-1 or clause 16 of the 3rd Ed of IEC 60601-1, respectively). Anybody connecting additional equipment to medical electrical equipment configures a medical system and is therefore responsible that the system complies with the requirements for medical electrical systems. Attention is drawn to the fact that local laws take priority over the above mentioned requirements. If in doubt, consult your local Medela customer service representative or your technical service department.

Please read and observe these warning and safety instructions before operation. These instructions for use must be kept with the device for later reference.

Please note that these Instructions for use are a general guide for the use of the product. Medical matters must be addressed by a physician.

Medela considers itself only responsible for the effect on BASIC SAFETY, reliability and performance of Thopaz+ if it is used in accordance with the Instructions for use.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

Subject to change.

WARNINGS

- “Pressure” generally implies “negative pressure”
- Consult this manual prior to use.
- Data transfer via USB is not allowed during therapy.
- After each use, the parts that have been in contact with aspirated secretions are to be cleaned and disinfected or disposed.
- Before cleaning the device, pull the plug out of the fixed mains socket.
- Thopaz⁺ was validated in combination with the accessories listed in Appendix A. For a correct and safe operation use Thopaz⁺ with these accessories only. Further information is supplied with the individual accessory.
- Do not use Thopaz⁺ in MRT (Magnetic Resonance Tomography).
- Do not dry Thopaz⁺ with microwaves.
- The device is not suitable for use while bathing, showering or in a hazardous explosive environment.
- The pressure range to be set must be determined by a physician in accordance with the age and weight of the patient.
- Do not use Thopaz⁺ if the drainage therapy indicates a pressure greater than max. pressure range of -10 kPa .
- Do not use Thopaz⁺ if the drainage therapy indicates a flow rate greater than max. flow capacity of 5 l/min .
- Do not use Thopaz⁺ if the drainage therapy indicates no pressure should be applied to the patient.
- Do not connect bilateral thoracic drains to one Thopaz⁺ unit. In such cases, the use of two Thopaz⁺ units is recommended.
- To ensure Thopaz⁺ can reduce the pressure, it is essential to wait 30 seconds between taking a drainage sample and unclamping the patient catheter (drain).
- For pediatric patients adapt pressure settings according to hospital guidelines.
- The use of Thopaz⁺ for any other indication than intended is neither desired nor allowed.
- For use only by medically trained and qualified persons who have been adequately trained in the use of Thopaz⁺.
- Wear gloves for all operations.
- The canister is replaced on the basis of a visual check or according to the instructions on the Thopaz⁺ display (warning signal).
- If a persistent air leak is shown, check the system is assembled correctly before taking further corrective actions. Ensure the system is air tight by clamping the catheter and observing the air leak decrease to zero.
- The catheter/connector interface is a location where clotting may occur. We recommend regular monitoring of this interface and appropriate removal procedure in case of occurrence.
- The Thopaz⁺ system is not designed for retransfusion.

- Do not touch the patient and the docking contacts simultaneously.
- This product contains lithium-ion batteries which bear risk of fire, explosion and burns. Do not disassemble, crash, heat above 100°C (212°F), incinerate or dispose of in fire.
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walki-talkie can affect the Thopaz⁺ and should be kept at least a distance 30 cm away from the equipment.
- HF (high-frequency) surgical equipment, radio networks or the like can influence the operation of the device and may not be operated in combination with the Thopaz⁺.

CAUTIONS

- No modification of this equipment is allowed.
- Before connecting Thopaz⁺ to the power supply, please verify that the supply voltage corresponds to that given on the device specification plate.
- Incorrect use of Thopaz⁺ can cause pain and injury to the patient.
- Do not use sterile accessories when the sterile packaging is damaged.
- Non sterile and reusable accessories must be cleaned and disinfected according to chapter 13 – General cleaning guidelines.
- For the initial startup do not switch Thopaz⁺ on before first charging the battery.
- The Thopaz⁺ 2l canisters are not intended for portable use (not by hand and not using the carrying strap).
- For safety reasons canister changes within the appropriate fluid level shall be recorded.
- Before diagnosing controll detected fluid value for plausability.
- For patients for whom a breakdown of the device can lead to a critical situation a replacement device must always be available.

Safety related tip

- Never place Thopaz⁺ with an attached 2l canister on the docking station without attaching the optional adapter docking station 079.0038 for Thopaz⁺ 2l canister first.
- The acoustic alarm “Leak in system – 301” is deactivated for approximately 6 minutes after Thopaz⁺ is switched on.
- Everytime a new patient is connected, it is recommended to carry out the functional check.
- Before the first application in accordance with the intended use, an inspection according to IEC 62353 (see Appendix B) is recommended to acquire the reference values.
- Medela recommends to use the biggest possible connector size compatible for the catheter.
- Use sample port of tubing to drain a sample. Therefore refer to chapter 11 – Taking a drainage sample.

Safety instructions

- Thopaz+ is a medical device that requires specific safety measures in regard to EMC. It must be installed and put into operation in accordance with the attached EMC information in chapter 19 – Technical documentation.
- In each of the following cases, Thopaz+ must not be used and it must be repaired by Customer Services:
 - If the power cord or the plug are damaged
 - If the device is not functioning according to routine check
 - If the device is damaged
 - If the device shows clear safety defects.
- Thopaz+ has no user serviceable parts inside. For safety reasons, it is required that Thopaz+ is repaired throughout its service life strictly and exclusively by Medela authorised service centres.
- Keep the power supply cord away from hot surfaces.
- The mains plug and the on/off switch must not come into contact with moisture. Never pull the mains plug out of the fixed mains socket by pulling on the power supply cord.
- Thopaz+ must stand upright during use.
- Separation from the mains is only assured through the disconnection of the mains adapter and fixed socket connection.
- Never use the device at high room temperatures, while bathing or showering, if you are very tired or in an environment where there is a risk of explosion.
- Never place Thopaz+ in water or other liquids.
- When attaching a new canister, verify that the canister size displayed on screen is the same size as the canister attached.
- When using single use, sterile products, please note that they are not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.
- Contact your local Medela customer service representative for assistance with product operations.
- The patient should be regularly monitored according to internal hospital guidelines.

These instructions for use must be kept for later reference.

2 Power supply and battery operation

WARNINGS

This product contains lithium-ion batteries which bear risk of fire, explosion and burns. Do not disassemble, crash, heat above 100°C (212°F), incinerate or dispose of in fire.

Thopaz+ can be operated from the mains power supply [REF](#) 087.0059, from the docking station [REF](#) 079.0037 or with the integrated lithium-ion battery. Before using the docking station, read first the docking station instructions [REF](#) 200.1554. The battery is charged during mains operation. The battery duration is dependent on the run-time of Thopaz+. This is influenced by the extent of parenchymal leakage and the set pressure. Thopaz+ does not run continuously but only switches on when the actual and nominal values differ.

During continuous operation, Medela guarantees a minimum run-time of 4 hours in battery operation after the device has been fully charged. In practical operation the real battery run time results in >10 hours.

If Thopaz+ is not used very often, the battery must be charged approximately every 6 months to ensure optimum functioning.

Backup battery

If an internal defect occurs (broken cable, battery defective) Thopaz+ turns off and the acoustic warning is sounded for at least 3 minutes (powered by a backup battery). Under these conditions Thopaz+ functions as a one way valve. Replace Thopaz+ immediately.

Charge battery

Plug the Thopaz+ power supply into a wall outlet and open the small protection cover to access the DC socket. Plug the coaxial power supply into the DC socket. Moving bars on the battery display will indicate charging is occurring [].

Battery low


Approximately 30 minutes before the battery is fully empty an acoustic signal sounds and the battery symbol [] starts to blink. The acoustic signal can be muted but the battery symbol continues to blink until the battery is empty. The set pressure is maintained but the battery should be recharged as soon as possible.

Press simultaneously to acknowledge an acoustic alarm and follow instructions on the display (see chapter 12 – Troubleshooting and alarms).

Battery empty

The battery symbol will blink with an acoustical sound for 10 minutes prior to the battery fully discharging, unless the pump is switched off earlier. The acoustic sound cannot be muted during the 10 minutes. If Thopaz+ is not connected to a power source the pump will switch off automatically after 10 minutes. The negative pressure is not maintained after the pump turns off.

1. Battery empty

2. Battery almost empty

3. Battery is being charged (moving bars)

4. Battery is fully charged and connected to the mains

* Symbol blinks

3 Description

Introduction

Thopaz+ is a high-quality digital chest drainage system. The compact system provides regulated negative pressure close to the patient's chest. It monitors the air leak and the fluid collection. Thopaz+ only applies vacuum required to maintain the negative pressure prescribed and set by the managing physician. The digital colour display provides objective data in real time as well as in historical graphs, which allows easy tracking of the therapy progress. It combines easy handling and cleaning with safety features to ensure optimal operation.

Intended use/indications

Thopaz+ is intended to be used for aspiration and removal of surgical fluids, tissue, gases, bodily fluids or infectious materials. Thopaz+ is indicated for all situations where chest drains are applied – especially for thoracic drainage in the pleural and mediastinal cavity in situations such as pneumothorax, after cardiac or thoracic surgery (post-operative), thorax injury, pleural effusion, pleural empyema or other related conditions. Thopaz+ is intended for use on patients in appropriate care settings.

Contraindications

There are no known contraindications for the Thopaz+ digital chest drainage system.

Intended user

Thopaz+ should only be operated by properly instructed users. These persons must not be hard of hearing or deaf and must have adequate visual faculty.

Intended patient population

Thopaz+ is intended to be used on patients only exhibiting conditions as described in the indications for use.

Important note

Compliance with proper surgical procedures and techniques is the responsibility of the physician. Each physician must evaluate the appropriateness of the treatment based on his own knowledge and experience.



4 Overview


WARNINGS

Do not touch the patient and the docking contacts simultaneously.

Main elements of Thopaz+

Display

How to activate a button



Press the buttons briefly to activate the desired action.

Exception: for the function switch to “Standby”, press the button for longer than 3 seconds.

Tubings

STERILE EO

Normal/small/large/single/double connector
Material: PVC (medical grade) sterile, twofold packaging
Length: 1.5 m

079.0025/26: Patient tube inner diameter: 5.5 mm
Connector outside diameter 1: 9.5 mm
Connector outside diameter 2: 16 mm

079.0021/22: Patient tube inner diameter: 5.5 mm
Connector outside diameter 1: 6.4 mm
Connector outside diameter 2: 16 mm

079.0023/24: Patient tube inner diameter: 5.5 mm
Connector outside diameter 1: 4.0 mm
Connector outside diameter 2: 11.0 mm

Thopaz tubing single,
large connector
REF 079.0025

Thopaz tubing double,
large connector
REF 079.0026

Thopaz tubing single
REF 079.0021

Thopaz tubing double
REF 079.0022

Thopaz tubing single,
small connector
REF 079.0023

Thopaz tubing double,
small connector
REF 079.0024

Safety related tip

Medela recommends to use the biggest possible connector size compatible for the catheter.

Canisters

STERILE EO


With/without solidifier 0.3l/0.8l/2l

Material: Polypropylene, sterile

Reading accuracy: +/- 2.5% (in upright position)

Solidifier: 0.3l = 12g / 0.8l = 28g / 2l = 2 x 30g

Simethicone capsule: Supports reliability of digital fluid measurement (0.3l = 720mg / 0.8l = 1440mg / 2l = 1440mg Simethicone)

Thopaz canister 0.8l, antifoaming

REF 079.0016

Thopaz canister 0.8l with solidifier, antifoaming

REF 079.0017

Thopaz canister 0.3l,
antifoaming

REF 079.0011

Thopaz canister 0.3l
with solidifier, antifoaming

REF 079.0012

Thopaz canister 2l,
antifoaming

REF 079.0018

Thopaz canister 2l
with solidifier, antifoaming

REF 079.0019

Safety related tip

Use sample port of tubing to drain a sample. Therefore refer to chapter 11 – Taking a drainage sample.

Safety related tip

The simethicone capsule requires 60 min of submersion for full effect.

5 Additional information

Definition of vacuum

By the application of medical aspiration devices, vacuum is normally given as the difference (in absolute figures) between absolute pressure and atmospheric pressure or as negative values in kPa. In this document, the indication of -10 kPa for example always refers to a pressure range in kPa below atmospheric ambient pressure (according to EN ISO 10079:1999).

Regulated pressure

The pressure is regularly checked and adjusted by Thopaz⁺. The set pressure corresponds to the pressure of the entire Thopaz⁺ system.

Safety pressure limit

The adjustable pressure range for Thopaz⁺ is -0,1 to -10 kPa (-100 mbar, -75 mmHg, -100 cm-H₂O). If pressure values greater than -7 kPa (-70 mbar, -53 mmHg, -70 cmH₂O) are set, the following warning appears in the display: too high pressure levels can cause pain and/or serious injuries to the patient. This warning must be acknowledged with "OK" before the pressure can be increased further.

Tube flushing

Flushing of the tubing occurs every 5 minutes or when Thopaz⁺ detects a siphon and prevents logging of the patient tubing.


Safety chamber in canister

For appropriate, safe operation Thopaz⁺ must be in an upright position. If Thopaz⁺ tips, the special construction of the safety chamber in the upper region of the canister protects the hydrophilic filter from premature clogging. The pressure is maintained. Important: immediately set Thopaz⁺ upright again. This function is disabled if the maximum filling capacity of the canister is exceeded.

Fluid measurement

Thopaz+ calculates the fluid amount from the fluid level in the canister and in regard to the attached canister size. The fluid level in the canister gets detected by capacitive coupling. For a valid fluid measurement Thopaz+ needs to stay within a tolerance of 10° vertical (upright) and the canister must contain a minimal fluid amount of:

Thopaz canister without solidifier	Thopaz canister with solidifier
30 ml in a 0.3l canister	40 ml in a 0.3l canister
50 ml in a 0.8l canister	70 ml in a 0.8l canister
200 ml in a 2.0l canister	

If these conditions are not fulfilled the fluid values appear in grey with a crossed out fluid symbol []. Additionally a Thopaz+ out of axis symbol [] or a Thopaz+ not enough fluid symbol [] appears on the action bar. No additional fluid can be measured in this status. The grey fluid values show the last measured fluid values.

This type of fluid measurement allows a measurement accuracy of +/- 15 ml in a 0.3l canister, +/- 40 ml in a 0.8l canister and +/- 100 ml in a 2.0l canister (corresponding to +/- 5% of the canister volume).

Thopaz+ displays the fluid values in a 0.3l canister in 1 ml steps, in a 0.8l canister in 5 ml steps and in a 2.0l canister in 10 ml steps. After a total drainage amount of 5000 ml and more the display will switch to liter (l). Thopaz+ will display fluid values in 0.01 l steps (equals 10 ml steps) thereupon.

Limitations

Thopaz+ calculates the fluid measurement according the fluid level in the respective canister. Movement of Thopaz+ can cause fluids to temporary cover the canister walls. A potential fluid film covering the canister wall can result in temporary higher fluid volume read out displayed on the screen. Also foam formation in the canister can result in higher fluid volume read out displayed on the screen. All Thopaz canisters (079.0011; 079.0012; 079.0016; 079.0017; 079.0018; 079.0019) are equipped with one or two simethicone capsules containing 720 mg simethicone each (see chapter 4 – Overview). Simethicone prevents foam formation and the coverage of fluid films on the canister wall.

Safety related tip

The simethicone capsule requires 60 min of submersion for full effect.

Safety related tip

For most accurate fluid measurement place or keep Thopaz+ on the docking station (079.0037).

Overflow protection / bacteria filter

A hydrophilic overflow protection/bacteria filter in the canister as well as in the tube sets (in the connector to Thopaz+) protects Thopaz+ from the entry of fluids and protects the environment from contamination.

Light sensor

Thopaz+ is equipped with a light sensor. It adjusts the display brightness to the surrounding light conditions. Thus, the patient does not get disturbed while sleeping by a bright shining display.

6 Installation

6.1 Check initial delivery

Check the delivery package of Thopaz+ for completeness and general condition.


Thopaz+

- REF** 079.1000 with Instructions for use 200.6841
- 079.1002 with Instructions for use 200.6842
- 079.1003 with Instructions for use 200.6843/101034256

Thopaz+/Liberty wall mains adapter/charger

- REF** 087.0059

Thopaz+ CD

- REF** 079.1007

Instructions for use

- REF** 200.6841 (for 079.1000)
- 200.6842 (for 079.1002)
- 200.6843/101034256 (for 079.1003)

6.2 Initial startup

WARNINGS



For use only by medically trained and qualified persons who have been adequately trained in the use of Thopaz+.

CAUTIONS

Before the first application in accordance with the intended use, an inspection according to IEC 62353 (see Appendix B) is recommended to acquire the reference values.

Important: Do not connect the patient to the system yet!

1. 
 - 1.1 Select plug.
 - 1.2 Attach plug to the wall mains adapter [REF 087.0059](#)
2. 
 - 2.1 Connect Thopaz+ to mains by mains adapter [REF 087.0059](#) or by docking station [REF 079.0037](#).
3. 
 - 3.1 Charge for approximately 3 hours or until the  symbol appears in the display.
 - 3.2 Press  to switch Thopaz+ on. A selftest starts.
4. 
 - 4.1 Select language.
 - 4.2 Confirm with "OK".
 - 4.3 Follow instructions on display.

5.

5.1 Press [] to switch Thopaz+ off.

6.

6.1 Disconnect Thopaz+ from the fixed mains socket by pulling on the plug housing. Do not pull on the cable or the bend protection.

7. Thopaz+ is ready for preparation for use.

7 Preparation for use/operating instructions

WARNINGS

- For use only by medically trained and qualified persons who have been adequately trained in the use of Thopaz+. Wear gloves for all operations. "Pressure" generally implies "negative pressure".
- Adapt pressure settings for pediatric patients.

CAUTIONS

- Do not use sterile accessories when the sterile packaging is damaged.
- Non sterile and reusable accessories must be cleaned and disinfected according to chapter 13 – General cleaning guidelines.
- Do not switch Thopaz+ on before first charging the battery.

7.1 Checks before use

- Check the Thopaz+ system before use for damage of the power cord or plug, obvious device damage or safety defects and proper functioning of the device.
- Check that the seal is correctly placed in the vacuum port of Thopaz+. If you do not know where the vacuum port is placed please refer to chapter 4 – Overview.
- Make sure that the rechargeable battery is recharged in case Thopaz+ is operated in battery mode.
- Check all accessories prior to use:
 - Canister for cracks, brittle and flawed spots.
 - Tubing for cracks, brittle areas and that connectors are firmly attached.
 Replace if necessary.


7.2 Connect tubing

1.

- 1.1 Open external packaging.
- 1.2 Keep patient connector in internal bag for hygienic reasons.

2.

- 2.1 Insert the smaller of the two connectors (connection to pump) horizontally in the direction of the arrow.

7.3 Snap in canister

1.

- 1.1 Unpack the canister.
- 1.2 Position the openings at the top and attach the bottom part of the canister to Thopaz+.

2.

- 2.1 Push the canister into Thopaz+ until you hear a click.

7.3.1 Using 2l canister

1.

- 1.1 Unpack the canister.
- 1.2 Position the openings at the top and attach the bottom part of the canister to Thopaz+.
- 1.3 Push the canister into Thopaz+ until you hear a click.

Safety related tip

Never place Thopaz+ with an attached 2l canister on the docking station without attaching the optional adapter docking station [REF 079.0038](#) for Thopaz+ canister 2l first.


CAUTIONS

The Thopaz+ 2l canisters are not intended for portable use (not by hand and not using the carrying strap).

7.4 Switch Thopaz+ on

Important: Do switch Thopaz+ on before connecting the system to the patient.

1. 
 - 1.1 Press [] to switch Thopaz+ on. A selftest starts.
 - 1.2 Check acoustic signal (beep).


If the selftest is not successful, follow the troubleshooting instructions on the display and refer to chapter 12 – Troubleshooting and alarms.

2. 
 - 2.1 Confirm whether a new patient is connected or not. The patient number is important for therapy time, history data and data transfer to the PC.

yes Thopaz+ will issue a new patient number. The therapy time and history data will start from 0 (recommended for a new patient).


Note: Once the patient number reaches 100 it automatically defaults back to 1.

no The patient number remains unchanged (recommended for the continued treatment of the same patient). The data history and the therapy time are retained.

3. 
 - 3.1 Check canister matching: detected canister size vs. attached canister. No match: dispose canister.

7.5 Carry out functional check

1.

1.1 Seal off the conical connector of the patient tubing with your thumb (through the packaging). The tubing must not be kinked otherwise the measuring tube will be closed.

2.

2.1 Press "OK"

Test passed:

Thopaz+ is ready for use.

Test failed:

Leak in the system:

- Check tubing connection.
- Check seal is correctly placed.
- Check pressure relief valve for tight connection.
- Replace canister and/or tubing if necessary.

3. Thopaz+ is in Standby and is ready for use with the factory settings. (To change factory settings please refer to chapter 10 – Changing settings)

4. Connect the Thopaz+ system to the patient catheter (drain) according to hospital guidelines. To start the therapy press "On".

Safety related tip

The acoustic alarm "Leak in system – 301" is deactivated for approximately 6 minutes after Thopaz+ is switched on.

Safety related tip

Everytime a new patient is connected, it is recommended to carry out the functional check.

7.6 Check therapy progress

7.6.1 Main display

The current size of the air leak (ml/min) is shown on the left side. The total amount of collected fluid (ml) during the therapy is shown on the right.

According to your preference the air leak history of 12h/ 100 ml/min and two customized fluid parameters can be displayed in addition. (To change factory settings please refer to chapter 10 – Changing settings).

- Air leaks between 0 and 1000 ml/min are displayed in 10 ml steps.
- Air leaks over 1000 ml/min are displayed in 100 ml steps.
- Rounding function
 - The digital display of air leak values is rounded to zero ml/min (0 ml/min) if the air leak measurement indicates an air leak < 5.0 ml/min. If the air leak measurement indicates an air leak \geq 5.0 ml/min and < 15 ml/min the digital display of air leak values is rounded to ten ml/min (10 ml/min)
- The fluid measurements are displayed as follow: In a 0.3 l canister in 1 ml steps, in a 0.8 l canister in 5 ml steps and in a 2.0 l in 10 ml steps. After a total drainage amount of 5000 ml and more the display will switch to liter (l). Thopaz+ will display fluid values in 0.01 l steps (equals 10 ml steps) thereupon.

7.6.2 Reset fluid display

For an individual fluid management an additional fluid parameter can be set to zero at any given time.


1.

1.1 Press [] for longer than 3 seconds.

2.

2.1 The timer on the left side is giving you the time when the parameter was set to zero (in this case 12 minutes ago). The fluid amount on the right side is giving you the drained fluid volume (in this case 50 ml got drained in the last 12 minutes). To reset the timer and fluid again press [] for longer than 3 seconds.

7.6.3 History


- Thopaz+ provides air leak history graphs and fluid history graphs.
- For scrolling through the graphs press "Next".
- To get back to your main display press "Home".

Air leak history

Press [] to get into the air leak history.

- Three different air leak history graphs (72 h/Autoscale, 24 h/1000 ml/min, 24 h/100 ml/min) plus the catheter check (see chapter 7.7 – Catheter check) are available.
- The air leak autoscale history graph adjusts the air leak scale to the largest measured air leak from the last 72 h to 1000/2000/3000/4000/5000 ml/min.



Fluid history

Press [] to get into the fluid history.

- Three different fluid history graphs are available.
72 h/Autoscale, 24 h/100 ml and 6 h/Autoscale.
- The fluid history graph shows the total collected fluid amount over the last 72 h.




7.7 Catheter check

Important: the catheter check works only with a patient catheter (drain) in the interpleural space.

Press [] to get into the air leak history. Scroll to page 4/4 by pressing "Next" [].

- The catheter check is only active if the air leak is 0ml/min.
- Follow instructions on screen.
- A pressure difference from minimal 3 cmH₂O between inspiration [] and expiration [] is displayed at an air leak of 0ml/min. This pressure difference is an indication for the respiration cycle of the patient and confirms that the patient catheter (drain) is not clogged.

▲ Patient catheter (drain) open

▼ Patient catheter (drain) clogged?

This number shows the measured maximal and minimal pressure of each breathing cycle. This pressure is indicated in the set pressure unit.

7.8 Switching Thopaz+ off

1.

- 1.1 Clamp patient tubing with the tube clamp.
- 1.2 Clamp patient catheter (drain).
- 1.3 Press “Standby” for longer than 3 seconds to switch off the pressure.

2.

- 2.1 Press [] to switch off Thopaz+.

3. Release, remove and seal the canister by the canister seal. Dispose of canister and patient tubing in accordance with internal hospital guidelines.
4. Clean Thopaz+ according to chapter 13 – General cleaning guidelines.

7.9 Transferring data to PC with ThopEasy+

With ThopEasy+ all stored data can be transferred to a PC as documentation, and for completion of the patient records. The data can be expanded with patient information, saved and printed. Operating instructions and the software can be found on the Thopaz+ CD [REF 079.1007](#) included in the scope of delivery.

For the connection between Thopaz+ and a PC, use the USB cable [REF 079.0034](#).

ThopEasy+ is compatible with the following Microsoft Operating Systems:

- Windows 7 32 Bit
- Windows 7 64 Bit
- Windows XP 32 Bit (SP3 or higher)
- Windows XP 64 Bit (SP3 or higher)
- Windows 8 32 Bit
- Windows 8 64 Bit
- Windows 10

WARNINGS

Data transfer via USB is not allowed during therapy.

8 Replace canister

WARNINGS

The canister is replaced on the basis of a visual check or according to the instructions on the Thopaz+ display (warning signal). Wear gloves for all operations.

CAUTIONS


Canister changes shall be recorded for safety reasons.

1. Prepare sterile canister (with/without solidifier).

2. Clamp off patient tubing with the tube clamp.

3.1 Switch Thopaz+ to Standby by pressing "Standby" for longer than 3 seconds.

4.1 Release the canister by pressing the release button.
4.2 Remove canister.

5.1 Unpack the canister.
5.2 Position the openings at the top and attach the bottom part of the canister to Thopaz+.

6.1 Push the canister into Thopaz+ until you hear a click.

Safety related tip

Using a 2l canister refer to chapter 7.3.1
– Using 2l canister.

7.  7.1 Check canister matching: Detected canister size vs. attached canister.
No match: Dispose canister.
- 7.2 Press "On". Pressure is built up and unclamp the tube clamp.
8.  8.1 Check credibility of air leak value.
9.  9.1 Seal the used canister with the canister seal.
10. Dispose of the used canister in accordance with internal hospital guidelines.

Canister with solidifier: activation of the solidification process

1.  1.1 Check if canister is sealed with the canister seal.
1.2 Press – to open the solidifier chamber.
2.  2.1 Shake – to activate the solidification process.
2.2 Dispose of the used canister in accordance with internal hospital guidelines.

9 Changing pressure during operation

WARNINGS

The negative pressure may only be changed by the physician or on medical orders.

Change pressure

Important: Thopaz+ is running.

Press simultaneously

- 1.1 Press [▶] and [◀] simultaneously.
- 1.2 Change pressure by pressing [⊗] or [⊗] and confirm with "OK" [◀].

9.1 Physiological pressure

For patients who are to be treated by gravity drainage (= water seal), the physiological pressure can be activated.

- 1.1 Press "Physio" [▶].
- 1.2 Confirm with "OK" [◀].

This mode corresponds to a pressure of -0.8kPa/
-6 mmHg/-8cmH₂O/-8 mbar.

10 Changing settings

10.1 Factory settings

Thopaz+ is delivered with the following factory settings. These settings can be changed and saved as the new standard.

Pressure:	-2.0 kPa – Choose to max. -10 kPa (-100 cmH ₂ O, -100 mbar, -80 mmHg).
Pressure unit:	kPa – Choose between kPa, cmH ₂ O, mmHg, mbar.
Measured pressure:	Hide – Choose between hide and show (only relevant in the air leak history). – This feature is added for advanced Thopaz+ users. A higher measured pressure than set pressure can be an indication that the patient is able to regulate the negative pressure in the inter-pleural space.
Fluid time 1:	OFF – Choose between 24 h / 12 h / 8 h / 4 h / OFF – Shows the total amount of collected fluid over the last X hours.
Fluid time 2:	OFF – Choose between 2 h / 1 h / 0.5 h (30 min) / 0.25 h (15 min) / OFF – Shows the total amount of collected fluid over the last X hours/minutes.
Fluid alarm:	OFF – Choose between 500 ml / 5min, 15min or 60min; 400 ml / 5min, 15min or 60min; 300 ml / 5min, 15min or 60min; 200 ml / 5min, 15min or 60min; 100 ml / 5min, 15min or 60min; 50 ml / 5min, 15min or 60min / OFF – The Fluid alarm 407 appears if Thopaz+ detected a higher fluid collection in the last 5 min, 15min or 60min than the desired parameter.
Air leak history:	Hide – Choose between hide and show. – This history can be displayed on the main display (see chapter 7.6.1 – Main display).
Language:	English – Choose your language.
Save as standard:	Choose if you want to save the new parameter as a standard (guidance see chapter 10.3 – Save a changed parameter as the new standard).

10.2 Change settings

WARNINGS

The settings may only be changed by the physician or on medical orders.

CAUTIONS

Adapt pressure settings for pediatric patients.


Important: Thopaz+ is on Standby.

1.

1.1 Press "Menu" [].

2.

2.1 Choose the desired parameter by pressing [].

2.2 Confirm with "OK" [] to change the desired parameter.


2.3 Change the desired parameter by pressing [] or [] and confirm with "OK" [].

2.4 Press "Home" [] to leave the settings.

2.5 The settings are changed for the current therapy.

10.3 Save a changed parameter as the new standard

1. Change all desired parameters as mentioned above.

2. Choose "Save as standard" and press "OK" []. "Saved" will appear for approx. 3 seconds.

3. The parameters are saved as the new standard. Whenever you switch Thopaz+ on and choose a new patient, Thopaz+ will be set with these settings.


4. Press "Home" [] to leave the settings.

11 Taking a drainage sample

WARNINGS

To ensure Thopaz+ can reduce the pressure, it is essential to wait 30 seconds between taking the sample from the sample port and unclamping the patient catheter (drain).


1. Before taking a drainage sample make sure that there is fluid in the patient tube.

2. Clamp the catheter (drain).

3.1 Switch Thopaz+ to Standby by pressing "Standby" for longer than 3 seconds.

4. Disinfect the sample port before taking the sample with CaviWipes or Mikrozid AF Wipes

5.1 Remove with syringe* air from the patient tubing. Repeat until fluid has gathered at the sample port.

* Syringe 17 G (1,4mm) or thinner.

6. With a syringe* take a sample from the patient tube.

7.1 Switch Thopaz+ on by pressing "On" – pressure is built up.

8. **Wait 30seconds!**

To ensure Thopaz+ can reduce the pressure, it is essential to wait 30 seconds between taking the sample and unclamping the patient catheter (drain).


9. Unclamp the patient catheter (drain).

12 Troubleshooting and alarms

Description

Thopaz+ distinguishes between warnings (yellow), alarms (red) and internal errors (red). If Thopaz+ detects any of these situations, an acoustic signal sounds and a description of the problem appears on the display. By pressing the two selection buttons simultaneously, the acoustic alarms are suppressed for 60 seconds. Exceptions: For the warning 306 "canister full" and the alarms 302 "tubing clogged", 313 "filter clogged" and 315 "Thopaz+ overheated" the acoustic alarms are suppressed for 5 minutes.

Example:

Visible on display					
	Error number	Problem description	Troubleshooting instructions	Remarks / Potential source of error	Pressure
Warning	306	Canister full	1. Clamp tubing 2. Switch to Standby 3. Replace canister 4. Continue with "On" 5. Unclamp tubing		Yes
	401	Battery low	Connect Thopaz+ to mains	Remaining time of battery approx. 30min.	Yes
	402	USB connection not permitted	Unplug USB cable	USB connection not permitted during operation or in the vicinity of the patient.	Yes
	405	Standby	Switch Thopaz+ on or off	In Standby for 5minutes.	No
	406	Temperature too high	Prepare backup device	Do not expose Thopaz+ to heat sources.	Yes
	409	Fluid level sensor out of order	– Fluid measurement unavailable – Confirm and continue therapy		Yes
	410	Canister detection	1. Clamp tubing 2. Switch to Standby 3. Replace canister 4. Continue with "On" 5. Unclamp tubing		Yes
Alarm	301	Leak in system	Check system for leaks	– Disconnection. – Missing seal at vacuum port.	No
	302	Tubing clogged	1. Clamp catheter 2. Clamp tubing 3. Switch to Standby 4. Replace tubing 5. Continue with "On" 6. Unclamp catheter	– Tubing is kinked or clogged for approx. 10min. – Thopaz+ detects clotting of Thopaz+ tubing, NOT of the patient catheter.	No
	305	Battery empty	Connect Thopaz+ to mains		No
	311	Selftest failed	Snap the canister out and in again	Do not switch Thopaz+ on if the patient is already connected.	No
	313	Filter clogged	1. Clamp tubing 2. Replace canister 3. Continue with "On" 4. Unclamp tubing	Filter in canister clogged.	No
	315	Thopaz+ overheated	Replace Thopaz+		No
	407	Fluid alarm	Patient is losing a large amount of fluid		Yes

Visible on display				
Error number	Problem description	Troubleshooting instructions	Remarks/Potential source of error	Pressure error
Internal error	Internal error XXX	1. Switch Thopaz+ off and on again 2. Inform Medela Customer Service	If error message repeats, inform Medela Customer Service.	

Safety related tip

The acoustic alarm “Leak in system – 301” is deactivated for approximately 6 minutes after Thopaz+ is switched on.

Safety related tip

Warning “canister full”/Alarm “filter clogged”.

The warning “canister full” is triggered when the canister filling level approximately reaches the maximal scale line indicated on the respective canister label. Change the canister according to the instructions of use, chapter 8 – Replace canister. If secretions line the walls of the canister, the warning “canister full” may be activated prematurely.

After muting the warning with the selection buttons, a symbol [✖] appears during operation.

After 5 minutes the warning reappears acoustically.

An overfilling of the respective canister volume will trigger the alarm “filter clogged”. The alarm is triggered due to the clogging of the hydrophilic filter upon contact with liquids or secretions.

Please refer to chapter “4 – Overview” for location of the hydrophilic, overflow protection/ bacteria filter in the canisters. Subsequent to a clogged filter and the alarm “filter clogged” the vacuum applied to the patient is interrupted and cannot be regulated or maintained.

An additional accumulation of liquids, secretions and air is released through the pressure relieve valve at the positive pressure of 0.2 - 0.5 kPa. This pressure release function is not a substitute for a free overflow function. Therefore the canister must be replaced once full. Please refer to chapter “4 – Overview” for location of the pressure relieve valve.

Change the canister according to the instructions for use, chapter “8 – Replace canister”.

13 General cleaning guidelines

WARNINGS

After each use, the parts that have been in contact with aspirated secretions are to be cleaned and disinfected or disposed.

General notes

- For detailed information refer to Medela Cleaning & Reprocessing Guide: 200.2391.
- This cleaning guideline relates to the cleaning of the entire Thopaz+ product ranges see Appendix A.
- Never place electrical devices in water or other liquids, specifically Thopaz+ and the docking station.
- Do not spray or pour liquid directly onto Thopaz+/docking station.
- Thopaz+ product range cannot be sterilized.
- Rinsing the Thopaz+ product range in a washing machine is not permitted.
- Immersion disinfection, thermal disinfection and ultrasound cleaning for the Thopaz+ product range are not permitted.
- Higher temperature than 45°C, 113°F can cause protein coagulation which can lead to problems further in the process.
- Clean all surfaces immediately after use to avoid residues from drying and to prevent growth of microorganisms.
- Do not use cleaning agent/detergent based on phenol.
- Do not use steel brushes or steel wool for cleaning.
- Store medical products dry and dust free.

Disposables

This symbol indicates a single use product. This product is not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.

1. Thorough cleaning of Thopaz[†]

Thorough cleaning can either be done at the point of use with sterile cold water only or in dedicated health care facilities with the additional use of enzymatic detergents according to the manufacturer's instructions for use.

- 1.1 Disconnect the power plug from the power socket.
- 1.2 Disinfect your hands and put on disposable gloves and proper protective gear.
- 1.3 Separate all individual parts. Dispose single use products in accordance with internal hospital guidelines.
- 1.4 Remove visible dirt with lint free nonwoven wipes wetted with sterile cold water (<40 °C, 104 °F). Clean around the vacuum port to ensure neither sterile water nor visible dirt enters the vacuum port.
- 1.5 Ensure that all surfaces are thoroughly wetted and keep moist for a minimum of 5 minutes then repeat this step with another wipe.
- 1.6 Dispose protective gloves and disinfect your hands.

2. Intermediate level disinfection

- 2.1 Disinfect your hands and put on disposable gloves and proper protective gear.
- 2.2 Use disinfecting wipes* according to the manufacturer's instructions for use. Clean around the vacuum port to ensure neither disinfecting agent nor visible dirt enters the vacuum port.
- 2.3 Let product work in for 5 minutes and wipe afterwards with another disinfecting wipe*.
- 2.4 Allow the surface to dry for a minimum of 5 minutes.
- 2.5 Dispose protective gloves and disinfect your hands.

*Recommended agent for intermediate level disinfection:

CaviWipes®

Metrex® Research

Address: 1717 West Collins Avenue

Orange, CA 92867, U.S.A.

Homepage: <http://www.metrex.com>

Phone: (800) 841 1428

Email: metrexcustcare@sybrondental.com

Mikrozid® AF Wipes

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Address: Robert-Koch Str. 2,

22851 Norderstedt, GERMANY

Homepage: <http://www.schuelke.com>

Phone: +49 (0) 40 521 00 0

Email: info@schuelke.com

3. Storage

- 3.1 After appropriate cleaning and disinfecting check the device before re-use.
 - 3.2 Store products dry and dust free under storing conditions (see chapter 17 – Technical specifications).
-

14 Warranty, maintenance and checks

Warranty

Medela AG warrants the device will be free from defects in materials and workmanship for a period of 2 years from the date of delivery ex works. Faulty material will be replaced free of charge during this period if not resulting from abuse or misapplication. This will not apply to parts subject to wear and tear in use. To ensure compliance with this warranty as well as optimum service from Medela products, we recommend the exclusive use of Medela accessories with our appliances.

In no event shall Medela AG be liable for claims which exceed the scope of warranty described including liability for consequential damages, etc. The right to the replacement of faulty parts will not be recognized by Medela if any work has been carried out on Thopaz+ by unauthorized persons. This warranty is subject to the appliance being returned to a Medela service centre.

Maintenance

Thopaz+ is a maintenance free device. Each time the device is switched on, a self-test is carried out which checks the internal functions of Thopaz+. A “beep” indicates that Thopaz+ successfully passed the self-test and that the device is switched on. Safe operation is ensured by the functional check including canister and tubing to be carried out before each use.

The safety instructions must be observed.

Routine check

Medela recommends subjecting Thopaz+ to a routine check once a year. The work to be carried out and a protocol template in accordance to EN/IEC 62353:2008 is given in Appendix B in these instructions.


Instructions for FPT test

Important: Patient must not be connected!
Thopaz+ must be disconnected from the mains (battery operation) and turned off.

Vacuum gauge
 REF 077.1456

REF 079.0021


REF 079.0016


(Tube e.g. $\varnothing 12 \times 7$ mm, not included in the scope of delivery)

1. Connect the tubing [REFI 079.0021](#) and the 0.8l canister [REFI 079.0016](#) to Thopaz+. For guidance to connect the system refer to chapter 7 – Preparation for use.
2. Connect vacuum gauge with tube to the patient connector.

3.

- 3.1 Press and hold [] and [] simultaneously. Press additionally [].
- 3.2 The safety test starts.
- 3.3 Follow the instructions on the display

Safety check

There are no prescribed or recommended safety checks to be carried out.

Rationales

- The engineering of the Medela Thopaz+ has been verified by independent test institutes for compliance with the EN / IEC 60601-1 standard. Copies of such certificates may be obtained on request. Because of the method of construction used, Medela does not expect electrical safety to be affected at any time during the life of the product – provided that Thopaz+ is repaired throughout its service life strictly and exclusively by Medela authorized service centers and that it is used properly in accordance with the intended use.
- For Thopaz+ as a device in protection class II (EN IEC 60601-1), there is no protective earth conductor; protective earth related tests (protective earth resistance, earth leakage current etc.) are not applicable.
- Thopaz+ enclosure is made entirely of insulating material. Tests of the touch current using common measuring instruments will therefore not reveal measurable values.
- The applied parts of Thopaz+ are connected to Thopaz+ via non-conductive vacuum hoses, canisters and tubings. Even when suctioning a conductive fluid (9 g/L NaCl; EN ISO 10079-1;1999) until the overflow protection device activates, measurements of the patient leakage current using common measuring instruments will not reveal measurable values.
- Thopaz+ systems do not have patient connections and functional earth conductors; tests related to those connections and currents are not applicable.

15 Disposal

Thopaz+ comprises metals and plastics and should be disposed of in accordance with the European directives 2011/65/EU and 2012/19/EU. Additional, local guidelines must also be observed. The electronic components must be disposed of separately, in accordance with the local regulations. This product contains lithium-ion batteries which bear risk of fire, explosion and burns, if disposed of improperly. Do not disassemble, crush, heat above 100 °C (212 °F), incinerate or dispose of in fire. Dispose of in accordance with your local requirements. Please take care that you dispose of Thopaz+ and its accessories in accordance with the hospital disposal guidelines.

User information for the disposal of electrical and electronic equipment

This symbol means that the electrical and electronic equipment must not be disposed as normal household refuse. A correct disposal of this device protects and prevents possible damage to the environment or human health. For more information about the disposal contact the manufacturer, your local caregiver or healthcare provider. This symbol is only valid in the European Union. Please respect the relevant state laws and rules in your country for the disposal of electrical and electronic equipment.

16 Accessories

WARNINGS

Thopaz+ was validated in combination with the accessories listed in Appendix A. For a correct and safe operation use Thopaz+ with these accessories only. Further information is supplied with the individual accessory.

17 Technical specifications

low vacuum,
-10kPa/-75mmHg/100cmHzO/mbar
Tolerance: +/-15%

Measured at 0m, atmospheric pressure:
1013.25 hPa
Please note: vacuum levels may vary
depending on location (meters above
sea level, atmospheric pressure and
temperature).

5l/min. (+/-10%)

1,0kg/2.2lbs


Power supply	Pump
Model: MSA-C2500IS12.0-30C-ZZ	12 VDC, 20W
IEC: 60601-1	
Input: 100-240 VAC, max 0.8 A, 50/60Hz	
Output: 12VDC, 2.5 A	


ISO 13485
CE (93/42/EEC), IIb

HxWxD
223x255x95 mm

Transport/Storage Conditions

Operating Conditions

IP33

18 Signs and symbols

This symbol indicates the compliance with the essential requirements of the Council Directive 93/42/ EEC of 14 June 1993 concerning medical devices.

This symbol indicates the compliance with additional USA and Canada safety requirements for medical electrical equipment.

This symbol indicates the legal specifications of the system.

This symbol indicates to follow instructions for use.

This symbol indicates to consult instructions for use.

This symbol indicates a CAUTION or WARNING associated with the device.

This symbol indicates safety related tip.

This symbol indicates the class of the system.

This symbol indicates a class II device.

IP33

This symbol indicates the protection against ingress of solid foreign objects and against harmful effects due to the ingress of water.

This symbol indicates a type CF applied part.

This symbol indicates a prescription device. U.S. Federal law restricts this device to sale by or on the order of a physician. (for US only).

This symbol indicates MR UNSAFE.

This symbol indicates that interferences may occur in the vicinity of equipment marked with this symbol.

This symbol indicates the manufacturer.

This symbol indicates the date of manufacture.

This symbol indicates that the device should not be used after the date shown.

This symbol indicates a single use device. Do not reuse the device.

This symbol indicates manufacturer's catalogue number.

This symbol indicates manufacturer's serial number.

This symbol indicates manufacturer's batch code.

This symbol indicates the device is sterilized using ethylene oxide.

This symbol indicates the temperature limitation for operation, transport and storage.

This symbol indicates the humidity limitation for operation, transport and storage.

This symbol indicates the atmospheric pressure limitation for operation, transport and storage.

This symbol indicates do not use the device if package is damaged.

unit

This symbol indicates a specified amount of quantity.

This symbol indicates to not dispose of the device together with unsorted municipal waste (for EU only).

This symbol indicates that the material is part of a recovery/recycling process.

This symbol indicates a carton package.

This symbol indicates to keep the device away from sunlight.

This symbol indicates to handle the fragile device with care.

This symbol indicates to keep the device dry.

This symbol indicates that audio is paused.

This symbol indicates the maximum vacuum level of the system.

This symbol indicates the flow levels of the system.

This symbol indicates the electrical specifications of the system.

This symbol indicates alternating current.

This symbol indicates direct current.

This symbol indicates the weight of the system.

This symbol indicates the dimensions (h x w x d) of the system.

pcs

This symbol indicates number of items.

This symbol indicates that device contains phthalates. Specific phthalates contained in the device are identified below the symbol.
DEHP: Di(2-ethylhexyl) phthalate

This symbol indicates air leak.

This symbol indicates fluid.

This symbol indicates Patient number.

This symbol indicates Thopaz+ is out of axis and fluid cannot be measured.

This symbol indicates to press the button more than 3 seconds.

This symbol indicates that battery is empty.

This symbol indicates that battery is almost empty.

This symbol indicates that battery is being charged (moving bars).

This symbol indicates that battery is fully charged and connected to the mains.

This symbol indicates fluid cannot be measured.

This symbol indicates not enough fluid for valid fluid measurement.

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19 Technical documentation

WARNINGS

Do not use other accessories than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased emissions or decreased immunity of the Thopaz+ device.

EMC

Thopaz+ is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and IEC 60601-1-2:2014 4th Edition according clause 7 and 8.9. Thopaz+ is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information.

Electromagnetic compatibility
 (EMC, IEC 60601-1-2:2007 3rd Edition and IEC 60601-1-2:2014 4th Edition, Table 1)

Electromagnetic emissions

Thopaz+ is intended for use in the electromagnetic environment specified below. The customer or the user of Thopaz+ should assure that they are used in such an environment.

Emission Tests	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	Thopaz+ uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	Thopaz+ is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

WARNINGS

Thopaz+ should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, Thopaz+ should be observed to verify normal operation in the configuration in which it will be used.

Electromagnetic compatibility
(EMC, IEC 60601-1-2:2007 3rd Edition, Table 2)

Electromagnetic immunity

Thopaz+ is intended for use in the electromagnetic environment specified below. The customer or the user of Thopaz+ should assure that they are used in such an environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 4 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV line(s) to earth not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of Thopaz+ requires continued operation during power mains interruptions, it is recommended that Thopaz+ be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels of a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Electromagnetic compatibility
(EMC, IEC 60601-1-2:2007 3rd Edition, Table 4)

Electromagnetic immunity

Thopaz+ is intended for use in the electromagnetic environment specified below. The customer or the user of Thopaz+ should assure that Thopaz+ is used in such an environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of Thopaz+, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 0.35 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 0.35 \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = 0.7 \sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed RF transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Thopaz+ are used exceeds the applicable RF compliance level above, Thopaz+ should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating Thopaz+.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

Recommended separation distance between portable and mobile RF communications equipment and Thopaz+

Thopaz+ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Thopaz+ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Thopaz+ as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m (meters)		
	150 kHz to 80 MHz $d = 0.35 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
0.01	0.04	0.04	0.07
0.1	0.11	0.11	0.22
1	0.35	0.35	0.70
10	1.11	1.11	2.21
100	3.50	3.50	7.00

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Appendix A

REF

- 079.1000 Thopaz+ digital cardiothoracic drainage system with Instruction for use 200.6841
- 079.1002 Thopaz+ digital cardiothoracic drainage system with Instruction for use 200.6842
- 079.1003 Thopaz+ digital cardiothoracic drainage system with Instruction for use 200.6843/101034256

Disposables

- 079.0011 Thopaz canister 0.3l, antifoaming*
- 079.0012 Thopaz canister 0.3l with solidifier, antifoaming*

- 079.0016 Thopaz canister 0.8l, antifoaming*
- 079.0017 Thopaz canister 0.8l with solidifier, antifoaming*

- 079.0018 Thopaz canister 2l, antifoaming*
- 079.0019 Thopaz canister 2l with solidifier, antifoaming*

- 079.0021 Thopaz tubing single*
- 079.0022 Thopaz tubing double*
- 079.0023 Thopaz tubing single, small connector*
- 079.0024 Thopaz tubing double, small connector*
- 079.0025 Thopaz tubing single, large connector*
- 079.0026 Thopaz tubing double, large connector*

Accessories

079.1006 Carrying strap Thopaz+

079.0037 Docking station

079.0038 Adapter docking station for Thopaz+ canister 2l

079.0034 USB cable

079.0039 Thopaz sealing cap*

079.0036 Holder with standard rail

079.0040 Universal holder with standard rail

077.1456 Vacuum gauge

087.0059 Thopaz+/Liberty wall mains adapter / charger

Documents

079.1007 Thopaz+ CD

200.6841 Thopaz+ Instr. for use (EN/DE/FR/IT/NL/SV/ES/PT)
 200.6842 Thopaz+ Instr. for use (EN/DA/FI/NO/CS/PL/SK/EL)
 200.6843 Thopaz+ Instr. for use (EN/JA/KO/RU/TR)
 101034256 Thopaz+ Instr. for use (EN/ZH)

Spare parts

077.0117	Screw covers	200.8057	Holder tool
077.0119	Feet	099.0347	Bed holder
099.0343	Seal orange	099.0361	Standard rail holder
077.0124	Cover for ports	099.0348	Release button

* **STERILE EO**

Appendix B

EN/IEC 62353: 2008

Hospital	Thopaz ⁺ serial number:
	Firmware version:
Location:	Run time (Settings 2/2): h

Visual checks	NO	YES
Cracks, brittle areas, discolouration on housing and display		
Covers for ports are present and can be clicked in		
Standard rail and bed holders are present and can be opened		
Seal is present, is clean and not brittle		
Release button can be pressed, canister can be clicked in		
Thopaz ⁺ is clean with no dirty areas		

Tests	OK	Not OK
FPT test Carry out test as described in Instructions for use. Thopaz ⁺ must be disconnected from the mains (battery operation). Required test equipment: vacuum gauge ~0-10 kPa Class 2.5% or better (E.G.077.1456)	<input type="checkbox"/>	<input type="checkbox"/>
Acoustic signal (beep) Audible?	<input type="checkbox"/>	<input type="checkbox"/>
Battery Connect Thopaz ⁺ to mains adapter or place it in docking station.	<input type="checkbox"/> Battery is being charged.	<input type="checkbox"/> Battery is not being charged.

Overall evaluation	Errors	Actions
No safety or functional technical errors were detected.	<input type="checkbox"/>	Thopaz ⁺ is ready for use.
Errors that are not safety relevant were detected. No direct risk.	<input type="checkbox"/> 1 _____ 2 _____ 3 _____	Thopaz ⁺ is ready for use, errors must be corrected by ___/___/___
Safety relevant errors were detected.	<input type="checkbox"/> 1 _____ 2 _____ 3 _____	Do not use Thopaz ⁺ until the errors are corrected. Send Thopaz ⁺ to Medela Service Centre.

Next test: ___/___/_____

Tester / Signature:	Place / Date:
---------------------	---------------

Notes

Notes

Medela AG
Lättichstrasse 4b
6340 Baar, Switzerland
www.medelahealthcare.com

International Sales

Medela AG
Lättichstrasse 4b
6340 Baar
Switzerland
Phone +41 41 562 51 51
Fax +41 41 562 51 00
customercare@medela.ch
www.medelahealthcare.com

Australia

Medela Australia Pty Ltd,
Medical Technology
3 Arco Lane
Heatherton, Vic 3202
Australia
Phone +61 3 9552 8600
Fax +61 3 9552 8699
contact@medela.com.au
www.medelahealthcare.com.au

Canada

Medela Canada Inc.
4160 Sladeview Crescent Unit # 8
Mississauga, Ontario
L5L 0A1
Canada
Phone +1 905 608 7272
Fax +1 905 608 8720
info@medela.ca
www.medelahealthcare.ca

France

Medela France SAS
14 rue de la Butte Corderie
91154 Etampes Cedex
France
Phone +33 1 69 16 10 30
Fax +33 1 69 16 10 32
info@medela.fr
www.medelahealthcare.fr

Germany

Medela Medizintechnik
GmbH & Co. Handel KG
Postfach 1148
85378 Eching
Germany
Phone +49 89 31 97 59-0
Fax +49 89 31 97 59 99
info@medela.de
www.medelahealthcare.de

India

Medela India Private Limited
602-603 | 6th Floor,
JMD Megapolis | Sohna Road,
Sector – 48 | Gurgaon 122 001
India
Phone +91 124 4832069
Fax +91 124 4832066
info@medela.in
www.medelahealthcare.in

Italy

Medela Italia Srl a socio unico
Via Turrini, 13–15 – Loc. Bargellino
40012 Calderara di Reno (BO)
Italy
Phone +39 051 72 76 88
Fax +39 051 72 76 89
info@medela.it
www.medelahealthcare.it

Japan

Medela K.K.
G-Place Sasazuka 13F
2-26-2, Sasazuka, Shibuya-ku,
Tokyo 151-0073
Japan
Phone +81 3 3373 3454
Fax +81 3 3373 3457
info@medela.jp
www.medelahealthcare.jp

Netherlands & Belgium

Medela Benelux BV
Uilenwaard 31
5236 WB's-Hertogenbosch
The Netherlands
Phone +31 73 690 40 40
Fax +31 73 690 40 44
info@medela.nl
info@medela.be
www.medelahealthcare.nl
www.medelahealthcare.be

Poland

Medela Polska Sp. z o.o.
ul. Wybrzeże Gdynskie
6D, 01-531 Warszawa
Poland
Phone +48 22 865 82 25
Fax +48 22 864 38 65
info@medela.pl
www.medelahealthcare.pl

Russia

OOO Medela
Vavilova str. 97, entrance 3
Office 16–17
117335 Moscow
Russia
Phone +7 495 980 6194
Fax +7 495 980 6194
info@medela.ru
www.medelahealthcare.ru

Spain & Portugal

Productos Medicinales Medela, S.L.
C/ Liacuna, 22 4º Planta
08005, Barcelona
Spain
Phone +34 93 320 59 69
Fax +34 93 320 55 31
info@medela.es
www.medelahealthcare.es
info@medela.pt
www.medelahealthcare.pt

Sweden

Medela Medical AB
Box 7266
187 14 Täby
Sweden
Phone +46 8 588 03 200
Fax +46 8 588 03 299
info@medela.se
www.medelahealthcare.se

Switzerland

Medela AG
Lättichstrasse 7
6340 Baar
Switzerland
Phone +41 848 633 352
Fax +41 41 562 14 42
contact@medela.ch
www.medelahealthcare.ch

United Kingdom

Medela UK Ltd.
Huntsman Drive
Northbank Industrial Park
Irlam, Manchester M44 5EG
United Kingdom
Phone +44 161 776 0400
Fax +44 161 776 0444
info@medela.co.uk
www.medelahealthcare.co.uk

USA

Medela LLC
1101 Corporate Drive
McHenry, IL 60050
USA
Phone +1 877 735 1626
Fax +1 815 307 8942
info-healthcare@medela.com
www.medelahealthcare.us